

## **Technical Datasheet**

**Analysis Name:** Free Amino Acids by UPLC

Method Number: LI-00.562

**Scope of Application:** This method applies to the determination of free amino acids

(except tryptophan), in protein hydrolysates, clinical nutrition, infant formulas, infant cereals and similar products, coffee/tea and related products, amino acid premixes, juice (except theanine and cystine), purees, pulps, and vegetable products

(e.g., carrots, sweet potatoes). It also applies to added

methionine in soy-based infant formulas.

**Description:** Aqueous extraction of free amino acids with or without

clarification using trichloroacetic acid (TCA). Determination of

amino acids by derivatization with 6-aminoquinolyl-N-

hydroxysuccinimidyl carbamate (AQC) using the AccQ•Tag Ultra Method (Waters Corporation, Milford, MA, USA). Derivatized amino acids are separated using reversed phase UHPLC with

UV detection at 260 nm.

Sample Weight 50 g

**Required:** 

**Analytical Platform: UPLC** 

**Special Information:** N/A

E ngacdublininfo@us.nestle.com

TDS-LI-00.562-5 9/5/2024



|                         |       |         | Limit of Quantification |                  |
|-------------------------|-------|---------|-------------------------|------------------|
| <b>Analyte Reported</b> | Alias | Unit    | Range                   | Reproducibility* |
| Histidine               | HIS   | mg/100g | 0.52 - 48.49            | 15%              |
| Asparagine              | ASN   | mg/100g | 0.44 - 41.25            | 15%              |
| Serine                  | SER   | mg/100g | 0.35 - 32.84            | 15%              |
| Glutamine               | GLN   | mg/100g | 0.49 -46.25             | 15%              |
| Arginine                | ARG   | mg/100g | 0.58 - 54.44            | 15%              |
| Glycine                 | GLY   | mg/100g | 0.25 - 23.46            | 15%              |
| Aspartic Acid           | ASP   | mg/100g | 0.44 -41.60             | 15%              |
| Glutamic Acid           | GLU   | mg/100g | 0.49 - 45.98            | 15%              |
| Threonine               | THR   | mg/100g | 0.40 -37.23             | 15%              |
| Alanine                 | ALA   | mg/100g | 0.30 - 27.85            | 15%              |
| G-Aminobutyric acid     | GABA  | mg/100g | 0.35 - 32.50            | 15%              |
| Theanine                | Thea  | mg/100g | 0.59 -55.00             | 15%              |
| Proline                 | PRO   | mg/100g | 0.38 - 35.98            | 15%              |
| Ornithine               | ORN   | mg/100g | 0.56 - 52.50            | 15%              |
| Cystine                 | CYS   | mg/100g | 0.40 -37.55             | 15%              |
| Lysine                  | LYS   | mg/100g | 0.49 - 45.69            | 15%              |
| Tyrosine                | TYR   | mg/100g | 0.60 - 56.63            | 15%              |
| Methionine              | MET   | mg/100g | 0.50 - 46.63            | 15%              |
| Valine                  | VAL   | mg/100g | 0.39 - 36.61            | 15%              |
| Isoleucine              | ILE   | mg/100g | 0.44 - 40.99            | 15%              |
| Leucine                 | LEU   | mg/100g | 0.44 - 40.99            | 15%              |
| Phenylalanine           | PHE   | mg/100g | 0.56 - 52.86            | 15%              |

<sup>\*</sup>Reproducibility: 17% for premixes, 20% for coffee components above 10mg/100g, and 40% for coffee components below 10mg/100g.

TDS-LI-00.562-5 9/5/2024