



NQAC

Nestlé Quality Assurance Center  
Dublin

# Technical Datasheet

**Analysis Name:** Determination of PAH by GC-MS/MS

**Method Number:** NQA-52.0005

**Scope of Application:** Oil, fat matrices, spices, herbs, teas, liquid milk, milk products as milk powder, skimmed milk powder and infant formula powder, cereals, and cocoa.  
Method is Not Applicable to gummies and white or black pepper extract.

**Description:** This method provides quantitative analysis of 4 priority polycyclic aromatic hydrocarbons (PAHs); benz[a]anthracene (BaA), benzo[a]pyrene (BaP), benzo[b]fluoranthene (BbF) and chrysene (CHR) in food by gas chromatography / tandem mass spectrometry (GC-MS/MS).

**Sample Weight Required:** 20g

**Method Reference:** EN 16619 Food analysis — Determination of benzo[a]pyrene, benzo[a]anthracene, chrysene and benzo[b]fluoranthene in foodstuffs by gas chromatography mass spectrometry (GC-MS).

**Analytical Platform:** GC-MS/MS

**Special Information:** PAHs are degraded by UV light. Protect PAHs solutions from light (wrap in foil if clear or translucent containers are used). Some plastic may have interactions with PAHs; use only polyethylene (PE) for sample containers if possible.



NQAC

Nestlé Quality Assurance Center  
Dublin

**Limits of Quantitation:**

| <b>Oil/Fat/spices/herbs/teas/<br/>and Liquid Milk</b> | <b>Milk powder/infant formula/cereal</b> | <b>Cocoa</b> |
|---|--|--------------|
| 1.0 ug/kg   | 0.1 ug/kg                                | 0.2 ug/kg    |

| <b>Analyte Reported</b> | <b>Alias</b> | <b>Reproducibility</b> |
|-------------------------|--------------|------------------------|
| benz[a]anthracene       | (BaA)        | ±25%                   |
| benzo[a]pyrene          | (BaP)        | ±25%                   |
| benzo[b]fluoranthene    | (BbF)        | ±25%                   |
| chrysene                | (CHR)        | ±25%                   |
| Sum of Quantified PAHs  | N/A          | N/A                    |